

# Light / Optical Microscopy (Metallography)

## Instrumentation and Sample Preparation

### Model:

- Nikon Epiphot Optical Microscope (Figure 1)
- The instrument uses light to study the microstructural properties of an alloy.

### Sample Preparation:

1. Sample is cut to approximately 15 mm.
2. The cut specimen is mounted in Bakelite to allow for easier sample preparation. (Figure 2)
3. Sample surface is finished using a series of sandpaper – from 120 grit up to 1200 grit. (Figure 3)
4. The sample is then polished using a 6, 3 then 1-micron cloth and diamond paste – to polish the surface to a mirror finish.
5. Finally – the sample is etched using mineral acids to reveal the microstructure, which can be examined using light microscopy.



Figure 1: Nikon Epiphot Optical Metallurgical Microscope





Figure 2: A Mounted and Polished Metallographic Sample




Figure 3: Silicon Carbide Paper (for finishing metallographic samples)




Figure 4: Microstructure of Sm<sub>2</sub>Co<sub>17</sub> (x10 magnification)

### The Microscope:

Light is generated from a light source (bulb), where it is passed through a series of lenses and focussed on the surface of the sample. The light is reflected off the sample surface and into an eyepiece and/or camera. This reveals the surface microstructure where micrographs can be taken for quality control (QC) purposes. (Figure 4)

The use of an inverted microscope gives greater freedom to move and manipulate the sample to thoroughly investigate the microstructure.

### Microstructural Assessment:

Typically, a sample is etched in a suitable etchant (e.g. nital, a mixture of nitric acid and alcohol). This reveals surface structure features e.g. grain boundaries (GB), non-metallic inclusions, precipitates etc. (Figure 5a/b, Figure 6)

Analysis can be done by using reference micrographs, graticules (eye-pieces) or using automated software.

## How does it work?

## What does it detect?

### Microstructural Properties:

- Grain size
- Metallic phases and secondary phases
- Grain distribution
- Segregation
- Defects
- Inclusions
- Porosity
- The volume fraction of phases
- Case depth and surface features

In addition to these features, it can be used as a QC tool to compare to previous melts of an alloy – as the microstructure can affect macroscopic, physical, mechanical and magnetic properties hence overall alloy performance under specific conditions.




Figure 5a

Figure 5b

Cross-section of Strip Cast NdFeB Flake micrographs (x10 and x20 magnification)




Figure 6: Microstructure of TiFeMn Powder Alloy